DRI OCT Triton[™] Series Swept Source Optical Coherence Tomography



NEW PixelSmart[™]



See. Discover. Explore.

The diagnostic power of Swept Source OCT Deep Range Imaging¹.

"

"Swept Source adds a new dimension to OCT. The TOPCON DRI Swept Source OCT is easy to use, provides unique clinical information, and has improved my practice. For the first time, we can in-vivo visualize not only the vitreo-retinal interface but also the cortical vitreous which is important at the time when more and more therapies are delivered via intra-vitreal injections. Deeper imaging brings choroidal thickness, helping guide my clinical decisions. Seeing more helps guide my therapy and allows me to treat more effectively. I find Swept Source OCT an essential tool to look for biomarkers of disease regression or progression."

Prof. P. E. Stanga

Manchester Royal Eye Hospital, Manchester Vision Regeneration (MVR) Lab at NIHR/ Welcome Trust Manchester CRF & University of Manchester

1) Fabio Lavinsky, Daniel Lavinsky. Novel perspectives on swept-source optical coherence tomography. Int J Retin Vitr (2016) 2:25

Welcome to the New Frontier in OCT Imaging.

The DRI OCT Triton combines the world's first² Swept Source OCT technology with multimodal fundus imaging.

Image Quality

Triton's Swept Source with its 100 kHz scanning speed and 1,050nm wavelength results in a clear and detailed images even for the deepest layers of the eye with short acquisition time. Visualize not only the retina and vitreous, but also the choroid and sclera¹.

Diagnostic capability

Seeing deeper makes it possible to have a better understanding of many ocular pathologies¹.With features such as OCT angiography, fundus autofluorescence and en face OCT, Triton empowers clinicians with multimodal imaging capability to help assess and preserve patient's eye health.


Practice efficiency

The Triton's automated functions, including single scan captures and SMARTTrack[™] system, are designed to optimize your practice workflow by simplifying data capture, analysis and follow-up.




See Deeper¹. See More.

Proliferative diabetic retinopathy




Courtesy: Prof. P. E. Stanga, Manchester Royal Eye Hospital, Manchester Vision Regeneration (MVR) Lab at N IHR/ Welcome Trust Manchester CRF & University of Manchester




Courtesy: Prof. P. E. Stanga, Manchester Royal Eye Hospital, Manchester Vision Regeneration (MVR) Lab at N IHR/ Welcome Trust Manchester CRF & University of Manchester

* FA photography and FAF photography can be performed using only DRI OCT Triton Plus.

Central serous chorioretinopathy



Courtesy: Prof. P. E. Stanga, Manchester Royal Eye Hospital, Manchester Vision Regeneration (MVR) Lab at N IHR/ Welcome Trust Manchester CRF & University of Manchester




Courtesy: Prof. P. E. Stanga, Manchester Royal Eye Hospital, Manchester Vision Regeneration (MVR) Lab at N IHR/ Welcome Trust Manchester CRF & University of Manchester

* FA photography and FAF photography can be performed using only DRI OCT Triton Plus.


See Deeper¹. See More.

Pathological myopia




Courtesy: Prof. P. E. Stanga, Manchester Royal Eye Hospital, Manchester Vision Regeneration (MVR) Lab at N IHR/ Welcome Trust Manchester CRF & University of Manchester

Macular pucker





Courtesy: Prof. P. E. Stanga, Manchester Royal Eye Hospital, Manchester Vision Regeneration (MVR) Lab at N IHR/ Welcome Trust Manchester CRF & University of Manchester

Image through cataract



a)





b)

a, b, c courtesy of Kazuya Yamagishi, MD (Hirakata Yamagishi Eye Clinic, Japan)

Discover What Lies Beneath

TOPCON's SS OCT Angio[™] combines OCT angiography with a Swept Source OCT. OCTARATM, a proprietary image processing algorithm, provides highly sensitive angiographic detection³, allowing for visualization of vascular structures even in the choroid and deeper retinal layers.

High-sensitivity Imaging and Deeper Intravascular Flow Visualization¹

Swept Source technology and OCTARA™ allow the deeper structures to be visualized with less depthdependent signal roll-off³. Additionally, the 1µm wavelength makes OCT Angiography imaging possible for patients with media opacities.

Rapid Scanning, Real Time Eye Tracking

At 100,000 A-Scans per second coupled with invisible^{*} scanning lines and the SMARTTrack[™] eye tracking system, 'the Triton quickly captures a dense data set and provides an en face OCT Angiography image of the retinal microvascular flow network'3.


* OCT Angiography scanning line may be faintly visible during capture to some people with certain conditions

Efficiency & Workflow Integration


Multimodal platform provides comparison of microvascular impairment with FA, FAF, OCT and true color fundus images in a single device*.


*DRI OCT Triton Plus

12x12mm 512 pixels



Graduate School of Medical Sciences, Hokkaido, Japan






icial 3.0x3.0

Courtesy: Akihiro Ishibazawa, MD, PhD, Asahikawa Medical University Graduate School of Medical Sciences, Hokkaido, Japan Retinal Research and Imaoing Austria

3) Magdy Moussa, Mahmoud Leila, Hagar Khalid. Imaging choroidal neovascular membrane using en face swept-source optical coherence tomography angiography Clinical Ophthalmology 2017:11 1859–1869


OCTARA[™]

OCTARA[™] is the image processing technology which extracts the signal changes derived from vascular flow using multiple OCT B-Scans acquired at the same position. It demonstrates high sensitivity for the detection of low blood flow in microvasculature³.




Multimodal Viewing

En face angiography images, B-scans and fundus photography can all be viewed on a single screen using IMAGEnet®6 and PinPoint[™] registration, so that area of interest can be assessed using multiple image modalities. Selected layers can easily be customized to enhance the clarity of specific pathological features.



Discover What Lies Beneath

Proliferative diabetic retinopathy




SS OCT Angio™ Montage

Courtesy: Akihiro Ishibazawa, MD, PhD.


Asahikawa Medical University Graduate School of Medical Sciences, Hokkaido, Japan

Before treatment




Courtesy: Akihiro Ishibazawa, MD, PhD. Asahikawa Medical University Graduate School of Medical Sciences, Hokkaido, Japan

After treatment



Courtesy: Akihiro Ishibazawa, MD, PhD. Asahikawa Medical University Graduate School of Medical Sciences, Hokkaido, Japan








Courtesy:Yuichiro Ogura, MD, Professor and Chairman, Department of Ophthalmology and Visual Science, Nagoya City University, Nagoya, Japan

Discover more possibilities: see beyond and deeper

OCTA metrics on Triton SS OCT Angio allows clinicians to objectively and quantitatively assess retinal vasculature, providing valuable insights into the patient's eye health.








Courtesy: Michael H. Chen, O.D.



Swept Source OCT Imaging

1,050nm wavelength


The longer wavelength light provides better tissue penetration, allowing visualization into the deepest layers of the eye¹.



Courtesy: Professor Jose Maria Ruiz Moreno MD, University of Albacete, Spain


Swept Source OCT technology; scanning speed of 100,000 A-scans/sec

The fast scanning speed of 100,000 A-scans/sec enables capture of clear B-scans⁴ by acquiring more A-scans within a given image acquisition time. This helps to reduce artifacts from involuntary eye movements such as saccades and blinks.



Invisible scan lines


The invisible 1,050nm wavelength light helps patients concentrate on the fixation target during the scan, reducing involuntary eye movement. It supports more efficient workflow in a practice by reducing the need to rescan.




Fabio Lavinsky, Daniel Lavinsky. Novel perspectives on swept-source optical coherence tomography. Int J Retin Vitr (2016) 2:25
Shoji Kishi. Impact of swept source optical coherence tomography on ophthalmology. Taiwan Journal of Ophthalmology 6 (2016) 58-68


En face OCT imaging

en face imaging allows for independent dissection of a depth range defined by two boundaries, selected from seven possible boundaries, by flattening the 3D data cube.





Courtesy: Prof. T. Nakazawa, MD,PhD, Tohoku University, Japan




Courtesy: Prof. T. Nakazawa, MD, PhD, Tohoku University, Japan

Vitreous visualization

Dynamic Focus[™]

Dynamic Focus[™] on Triton allows for acquisition of images with near uniform focus and image quality throughout the entire depth of the image, which helps to enhance the typically weaker signal from the vitreous.

> Clear image in all area



EVV (Enhanced Vitreous Visualization™)

EVV helps clinicians assess vitreous and vitreoretinal interface abnormalities¹. Contrast can be quickly adjusted to the needs of the physician, depending on the area of greatest interest.





Discover from Cornea to Choroid



Courtesy: Michael H. Chen, O.D.

Anterior segment imaging


Triton's (optional) anterior segment imaging capabilities allow for visualization of the cornea, anterior chamber angle, iris and sclera⁵.

Image samples


OCT image B-scan length 16mm



Anterior segment in Radial scan




Anterior segment in 3D scan



7 boundary segmentation/5 layer thickness maps/caliper function

Retinal tissue layers are automatically segmented by the Topcon Advanced Boundary Software (TABS[™]), enabling quantification of retinal thickness and sub layers⁶⁷.




Choroidal thickness maps

Understanding choroidal structural changes in relation to various disease and disease management, has been of long-standing interest to clinicians⁸. Triton provides clear visualization of the choroid, and generates choroidal thickness maps to aid in the assessment of choroidal structure.

Retina	between the ILM-OS/RPE boundaries
RNFL	between the ILM-RNFL/GCL boundaries
GCL+	between the RNFL/GCL-IPL/INL boundaries
GCL++	between the ILM-IPL/INL boundaries
CSI	between the BM-CSI boundaries or ILM-CSI boundaries

Hood Report (for Glaucoma)

Retinal Thickness/RNFL/GCL and Optic Nerve Metrics in just one scan. This report streamlines the decision-making process through the correlation of structure (GCC/RNFL) with function (overlay of visual field test locations)⁶.




- 6) Zhichao Wu, Denis S. D. Weng, Rashmi Rajshekhar, Abinaya Thenappan, Robert Ritch, Donald C. Hood. Evaluation of a Qualitative Approach for Detecting Glaucomatous Progression Using Wide-Field Optical Coherence Tomography Scans. Trans Vis Sci Tech. 2018;7(3):5.
 7) Beatriz Abadia, Ines Suñen, Pilar Calvo, Francisco Bartol, Guayente Verdes, Antonio Ferreras. Choroidal thickness measured using swept-source optical coherence
- Beatriz Abadia, Ines Suñen, Pilar Calvo, Francisco Bartol, Guayente Verdes, Antonio Ferreras. Choroidal thickness measured using swept-source optical coherence tomography is reduced in patients with type 2 diabetes. PLoS ONE 13(2): e0191977.
- Sushmitha Rao Uppugunduri, Mohammed Abdul Rasheed, Ashutosh Richhariya et al. Automated quantification of Haller's layer in choroid using swept-source optical coherence tomography. PLoS ONE 13(3):e0193324

DRI Meets Multimodal Fundus Imaging: See the Whole Picture

Swept Source OCT incorporates multimodal fundus imaging

DRI OCT Triton can acquire the OCT and fundus image in a single capture. PinPoint[™] Registration identifies the location of the B-scan on the fundus image. Comparison between the B-scan and fundus image can support clinical efficiency during diagnosis.



Courtesy: Jay M. Haynie, O.D.

OCT + Color fundus



Courtesy: Jay M. Haynie, O.D.

OCT + FAF

True color* Fundus images

The DRI OCT Triton offers a true color, non-mydriatic fundus image. Fluorescein Angiography (FA) and Fundus Autofluorescence (FAF) are available** to enhance the diagnostic capability of Triton Plus. The all-in-one device supports efficient workflow in practice.

*Color fundus image with white light, with 24-bit color. **DRI OCT Triton Plus :

OCT /Anterior OCT (Option)/ OCT Angiography (Option) /Color / Red-Free / FA / FAF

DRI OCT Triton : OCT /Anterior OCT (Option)/ OCT Angiography (Option) /Color/ Red-Free

Sensor Information

Color Imege Sensor :	5MP
FA Image Sensor :	4MP
FAF Image Sensor :	4MP




Stereo photography

Three dimensional visualization of color fundus images can be achieved by acquiring images in stereo photography mode. Triton's on-screen acquisition guidance supports quick and easy operation with auto alignment for capturing stereo pairs.



Panoramic wide field photography

In addition to macular and disc imaging, the Triton provides wide coverage of the retina. A panoramic graphic can be created from multiple fundus or OCT Angiography images.






Smarter Tracking. Smarter Workflow.

SMARTTrack[™]


The SMARTTrack[™] system enhances the tracking and follow-up ability of the Triton with a variety of functions designed to enhance its user-friendliness:

- Fundus-Guided Acquisition (FGA)
- Follow-up Function
- Tracking Photography




Fundus Guided Acquisition (FGA)

OCT scan location can be easily set by selecting the scan area on the fundus image. With FGA, the operator can choose to capture or import a fundus image, select the scan location and automatically acquire a B-scan at that location.



Follow-up Function

This function allows you to retrieve and re-analyze the same location at follow-up, for comparison of past and current images. All an operator needs to do is simply select the past data and Triton automatically captures the same area.





Motion Correction/Compensation/ Rescanning Function

Motion Correction


Corrects the Z direction movement


Compensation Function

Tracks the eye and then compensates for the X direction movement.

Rescanning Function

Parts of the scanning area may be missed due to Y direction eye movement. In such a case, the rescanning function automatically activates to rescan the missing area.





Before compensation

After compensation

Alignment navigation

When an operator wishes to acquire an image, the Triton's monitor guides the operator to reduce potential errors and make operation simple.

- Auto focus and auto shoot, in color/FAF mode
- Auto focus, auto-Z and Z-lock function, in OCT mode







too close

too far

Small pupil solution

Live Fundus View

The fast scanning speed allows the Triton to create a live en face fundus image, an ideal tool for precisely visualizing the scan position. This makes the disc, retinal vessels and scanning position easy to see, even in patients with small pupils.

OCT capture mode without retinal photography

Triton can also capture a 3D scan, with or without color fundus photography, to avoid a miotic response and better meet the needs of patients with the smallest pupils.



Reinventing Swept Source Imaging:

Triton with PixelSmart[™] is the next stage in Swept Source OCT imaging


NEW! PixelSmart[™]

Triton with PixelSmart The next level in Swept Source OCT Imaging

PixelSmart is Topcon's new image processing algorithm which reduces speckle noise to improve contrast, producing a better image quality.

PixelSmart is available for all Triton 3D scans, existing and new:

- 3D Wide
- 3D Macula
- 3D Disc
- Combination scans





Specifications

Observation & Photography of Fundus Image			
Photography Type	Color, FA*, FAF*, Red-free**		
Picture Angle	45°		
Picture Angle	43 Equivalent 30° (Digital Zoom)		
Operating Distance	34.8mm		
Photographable Diameter of Pupil	Normal: \$4.0mm or more		
	Small pupil mode diameter: \$\$.3mm or more		
Observation & Photography of Fundus			
Scanning Range (on fundus)	Horizontal Within 3 to 12mm		
Scalling Range (on fundus)	Vertical Within 3 to 12mm		
Scan Pattern	3D scan		
Scan Pattern			
Coop Croood	Linear scan (Line-scan/Cross-scan/Radial-scan/Raster-scan)		
Scan Speed Lateral Resolution	100,000 A-Scans per second		
	20µm Disital 2 Cum		
In-depth Resolution	Digital: 2.6µm		
Dhotographable Dispector of Dupil	Optical function: 8µm		
Photographable Diameter of Pupil	¢2.5mm or more		
Observation & Photography of Fundus			
Fixation target	Internal fixation target :		
	Dot matrix type organic EL		
	The display position can be changed and adjusted.		
	The displaying format can be changed.		
	Peripheral fixation target :		
	This is displayed according to the internal fixation target		
	indicated position.		
	External fixation target		
Observation & photography of anterior			
Photography type	IR		
Operating distance	17mm		
Observation & photography of anterior			
Operating distance	17mm		
Scan range (on cornea)	Horizontal Within 3 to 16mm		
	Vertical Within 3 to 16mm		
Scan pattern	3D scan		
	Linear scan (Line-scan/Radial-scan)		
Scan speed	100,000 A-Scans per second		
Fixation target	Internal fixation target		
	External fixation target		
Electric Rating			
Power Source	Voltage: 100-240V		
	Frequency: 50-60Hz		
Power input	250VA		
Dimensions / Weight			
Dimensions	320-359 mm(W) X 523-554 mm(D) X 560-590 mm(H)		
Weight	21.8kg (DRI OCT Triton)		
	23.8kg(DRI OCT Triton Plus)		
* EA whether we what and EAE whether we what are the weather would be the	with the DDL OCT Tritee Dlue		

 * $\,$ $\,$ FA photography and FAF photography can be performed only with the DRI OCT Triton Plus.

** In this digital red-free photography, the color image is processed and is displayed as a pseudo red-free photographed image.

*** Observation & photography of anterior segment can be performed only when the anterior segment attachment is used.







TOPCON CORPORATION

75-1 Hasunuma-cho, Itabashi-ku, Tokyo 174-8580, JAPAN, Phone: +81-(0)3-3558-2522/2502 Fax: +81-(0)3-3965-6898 www.topcon.co.jp

TOPCON MEDICAL SYSTEMS, INC. 111 Bauer Drive, Oakland, NJ 07436, U.S.A. Phone: +1-201-599-5100 Fax: +1-201-599-5250 www.topconhealthcare.com

TOPCON CANADA INC. 110 Provencher Avenue, Bolsbriand, QC J7G INI CANADA Phone: +1-450-430-7771 Fax: +1-450-430-6457 www.topcon.ca

TOPCON HEALTHCARE SOLUTIONS, INC. 111 Bauer Drive, Oakland, NJ 07436, U.S.A U.S.A. Phone: 1-201-599-5100 Fax: 1-201-599-5250 www.topconhealthcare.com

TOPCON HEALTHCARE SOLUTIONS ASIA PACIFIC PTE. LTD. Jalan Kilang Timor, #09-01 Pacific Tech Centre, SINGAPORE 159303 Phone: +65-68720606 E-mail: medical_sales@topcon.com.sg www.topcon.com.sg

/ww.topcon.com.sg TOPCON HEALTHCARE SOLUTIONS EMEA OY HQ & Product Development Saaristonkatu 23 90100 Oulu, Finland Phone: +358 20 734 8190

www.topconhealth.eu

TOPCON HEALTHCARE SOLUTIONS TOPCON HEALTHCARE SOLUTIO AUSTRALIA PTY LTD 14 Park Way, Mawson Lakes, South Australia, 5095, Australia Phone: +61-8-8203-3306

Email: au.info@topcon.com www.topconhealth.com.au EC REP TOPCON EUROPE MEDICAL B.V. Essebaan 11, 2908 LJ Capelle a/d IJssel, P.O.Box 145, 2900 AC Capelle a/d IJssel, THE NETHERLANDS THE NETHERLANDS Phone: +31 -(0)10-4585077 Fax: +31 -(0)10-4585045 E-mail: medical@topcon.com www.topconhealthcare.eu

www.topconnealthcare.eu **TOPCON ITALY** Viale dell'Industria 60, 20037 Paderno Dugnano, (Miano), ITALY Phone: +39-02-91081091 E-mail: Info@topcon.it.

www.topcon-medical.it www.topcon-medical.it **TOPCON DANMARK** Praestemarksvej 25, 4000 Roskilde, DANMARK Phone: +45-46-327500 Fax: +45-46-327555 E-mail: Info@topcon.dk www.topcon-medical.dk

TOPCON IRELAND MEDICAL Unit 292, Block G, Blanchardsto Corporate Park 2 Ballycoolin Dublin 15, D15 DX58, IRELAND Phone: +353-12233280 dstown. E-mail: medical.ie@topcon.com nedical.ie www.topcoi

www.topcon-medical.ie **TOPCON DEUTSCHLAND MEDICAL G.M.B.H.** Hanns-Martin-Schleyer Strasse 41, D-47877 Willich, GERMANY Phone: +49-(0)2154-8850 Fax: +49-(0)2154-88517 E-mail: info@topcon-medical.de www.topcon-medical.de

TOPCON SCANDINAVIA A.B. Neongatan 2, P.O.Box 25, 43151 Mölndal, SWEDEN Phone: +46-(0)31-7109200 Fax: +46-(0)31-7109249 E-mail: medical@topcon.se

www.topcon-medical.se www.topcon-medical.se **TOPCON ESPAÑA S.A.** Frederic Mompou, 4, 08960 Sant Just Desvern Barcelona, SPAIN Phone: +34-93-4734057 Fax: +34-93-4733932 E-mail: medica@topcon.es www.topcon-medical.es

TOPCON (GREAT BRITAIN) MEDICAL LTD. Topcon House, Kennet Side, Bone Lane, Newbury, Berkshire, RG14 5PX, UK Phone: 444-1635-551120 Fax: +444-1635-551170 E-mail: medical@topcon.co.uk www.topcon-medical.co.uk TOPCON FRANCE MEDICAL S.A.S.

1 rue des Vergers, Parc Swen, Bâtiment 2, 69760 Limonest, FRANCE Phone: +33-4-37581940 Fax: +33-4-72238660 Email: topconfrance@topcon.com www.topcon-medical.fr

TOPCON POLSKA SP. Z. O. O. ul. Warszawska 23, 42-470 Siewierz, POLAND Phone: +48-(0)32-6705045 Fax: +48-(0)32-6705045 Fax: +48-(0)32-6713405 Email: info@topcon-polska.pl www.topcon-medical.pl

TOPCON SINGAPORE MEDICAL PTE. LTD. 1006 Pasir Panjang Road, #05-05, Interlocal Centre, SINGAPORE 118523 Phone:+65-68720606 Fax:+65-67736150 E-mail: med.sales.sg@topcon.com www.topcon.com.sg

TOPCON INSTRUMENTS (MALAYSIA) SDN.

BHD No.6, Jalan Pensyarah U1/28, Hicom No.6, Jalah Pensyarah O/28, mcom Glenmarie Industrial Park, 40150 Shah Alam, Selangor, MALAYSIA Phone: +603 5023 2688 Email : mys_medical_sales@topcon.com

TOPCON INSTRUMENTS (THAILAND) CO., LTD. OPCON INSTRUMENTS (THAILANI 77/162 Sinnsathorn Tower, 37th F1 Krungthonburi Rd., Klongtonsai, J Bangkok 10600, THAILAND Phone: 662-4401152 Fax: 662-4401158 Email : tha_medical@topcon.com ai. Klongsarn

MEHRA EYETECH PRIVATE LIMITED

8018 Wing, Lotus Corporate Park, Grah Firth Steel Compound Goregaon (East) Mumbai 400063 Maharashtra, INDIA Phone: +91-22-61285455 www.mehraeyetech.in Graham

TOPCON (BEIJING) MEDICAL TECHNOLOGY CO., LTD. Room 2808, Tower C, JinChang An Building, No.82, Middle Section of East 4th Ring Road, Chaoyang District, Beijing 100124, P.R. CHINA Phone: +86-10-8794-5176

IMPORTANT In order to obtain the best results with this instrument, please be sure to review all user instructions prior to operation. Not available for sale in all countries. Please check with your local distributor for availability in your country.



