## Topcon Slit Lamp Series





# The Combined Result of **Optics and Technology** SL-D series Digital Slit Lamps

The blending of quality, versatility and ease of use make digital image capture smooth and effortless. The SL-D series meets the varied demands of eye care specialists.



### Optical Quality

With more than 85 years of experience in lens manufacturing, Topcon produces all the optics to equip its instruments, ensuring excellent lens quality and a bright, clear view.



### Ergonomics

Topcon Slit Lamps are ergonomically designed for comfortable use. The SL "D" Series employ converging binocular tubes with a parallel magnification system, for comfortable viewing.



### Versatility

A variety of models in different configurations offering primary clinical observation through to sophisticated image capture. The Topcon series will satisfy any need and budget.



### LED Illumination

As a leader in ophthalmic instrumentation, Topcon offers LED as the recommended illumination in many of its Slit Lamp models. The lifetime of LED is substantially longer than conventional halogen bulbs. LEDs don't need frequent replacement, making Topcon Slit Lamps environmentally sound.


### **Filter System**



### Topcon Blue-Free Filter™ System¹

Observation can be performed using the Blue-Free Filter<sup>™</sup> System, by using the incorporated exciter filter with the optional barrier filter<sup>2</sup>. Tear dynamics and subtle ocular surface damage can be more clearly viewed, with greater contrast and fewer artefacts.<sup>4</sup>

- Optional accessory
   Shizuka Koh et al., "Diagnosing Dry Eye Using a Blue-free Barrier Filter", American journal of ophthalmology, 2013 September
   Images provided by Toru Noda, M.D, Department of Ophthalmology, Tokyo Medical Center
   Arita R, Morishige N, Fujii T, et al., "Tear Interferometric Patterns Reflect Clinical Tear Dynamics in Dry Eye Patients", Investigative Ophthalmology and Visual Science, July 2016



Blue-Free Filter<sup>3</sup>

Cobalt Blue Filter<sup>3</sup>



### **Amber Filter**

The SL-D Series incorporates an amber filter, which can be easily selected.








Non-Amber

### High quality digital images with the DC-4 Digital Photo Attachment\*

The DC-4 Digital Camera can be easily mounted on any of the "D" series Topcon Slit Lamps, converting them into an ophthalmic image recording system The DC-4 captures color and monochrome images of the eye. When used in conjuction with the BG-5 background illuminator, the infrared capabilities allow viewing and documenting of the Meibomian glands.









### **Features**



**Smart Capture Function** five images captured in rapid succession



**Live View** see exactly what you are capturing on screen



Still Image with Auto **Exposure (AE) Function** always the correct color and brightness balance



**Video Capability** video clips are easily obtained and displayed




**Meibomian Gland Status** view and document Meibomian glands condition with infrared viewing

### **Meibomian Gland Observation**


The Meibomian glands can be examined and their status documented using the optional BG-5 illumination system in selected models.\*1\*2

\*1 The selected configuration is Digital Slit Lamp SL-D701 + Digital Camera DC-4 + background illuminator BG-5.
\*2 Arita Reiko et., "Development of Definitive and Reliable Grading Scales for Meibomian Gland Dysfunction", American journal of ophthalmology September 2016;169:125-137




### **Image Examples**

The DC-4 digital sensor provides natural color and sharp, clear images.



### Connectivity When used in combination with the DC-4 Digital attachment, all SL "D" Series Slit Lamps can capture and transfer images to the IMAGEnet® 6 ophthalmic data system. IMAGEnet® 6 is Topcon's proprietary data management software capable of capture, storage, display and enhancement of digital images obtained with Topcon devices.

### **Topcon IMAGEnet® 6 Data Management System**



The balanced offering of LED and halogen illumination models plus a selection of optional accessories, provides versatility and convenience.

### **Slit Lamp - halogen type**



### SL-D2/D4

Digital Slit Lamp

**SL-D2:** A basic Slit Lamp with three magnifications, halogen illumination and digitally ready for use with the DC-4 attachment.

**SL-D4:** Five magnifications, halogen illumination and digitally ready for use with the DC-4 attachment.



### SL-D301

Digital Slit Lamp

Basic three magnifications Slit Lamp with halogen tower illumination system and digitally ready for use with the DC-4 attachment.



### **SL-D701**

Digital Slit Lamp

Halogen illumination, five magnifications and digitally ready for use with the DC-4 attachment. Accepts the BG-5 background illuminator with infrared capability for meibography capture and the Blue-Free™ Filter System for fluorescein observation.

### Slit Lamp - LED type



### SL-2G

LED Slit Lamp

Basic, economical and environmentally friendly with three magnifications and LED illumination.



### SL-D4 (Type: LED)

Digital Slit Lamp

Five magnifications, LED illumination and digitally ready for use with the DC-4 attachment. Includes exciter filter for use with optional Blue-Free  $^{\text{TM}}$  Filter observation system.




### **SL-D701**


Digital Slit Lamp

LED illumination, five magnifications and digitally ready for use with the DC-4 attachment. Accepts the BG-5 background illuminator with infrared capability for meibography capture and the Blue-Free™ Filter System for fluorescein observation.


### System Configuration: SL-D2/D4 (Type: LED )




### **System Configuration: SL-D301**



### **System Configuration: SL-D701**





 $<sup>^{\</sup>star}$  Not available in all countries. Please check with the distributor in your country for availability.

Microscope	SL-2G	SL-D2	SL-D4/ SL-D4(Type:LED)	SL-D301	SL-D701
Туре	Galilean type	Galilean type		Galilean type	Galilean type
Magnification Selection	3 steps	3 steps	5 steps	3 steps	5 steps
	by drum rotation	by drum rotation	by drum rotation	by drum rotation	by drum rotation
yepiece	12.5x	12.5x	12.5x	12.5x	12.5x
Aagnification Ratio	10x(22.5mm)	10x(22.5mm)	6x(35.1mm)	10x(22.5mm)	6x(35.1mm)
(Field of View)	16x(14.1mm)	16x(14.1mm)	10x(22.5mm)	16x(14.1mm)	10x(22.5mm)
	25x(8.8mm)	25x(8.8mm)	16x(14.1mm)	25x(8.8mm)	16x(14.1mm)
	,		25x(8.8mm)		25x(8.8mm)
			40x(5.6mm)		40x(5.6mm)
Slit Illumination Slit Width	0 -14mm, can be altered gradually (14mm=circle)	0 -14mm can be	altered gradually (14mm=circle)	0 - 9mm, can be altered gradually (9mm=circle)	0 -14mm, can be altered gradually (14mm=circle)
Slit Length	1-14mm, can be altered gradually (14mm=circle)		oltered gradually (14mm=circle)	1-8mm, can be altered gradually	1-14mm, can be altered gradually (14mm=circle)
Filters	Blue. Red-Free		e. Red-Free. Amber	Blue, Red-Free	Blue filter, Red-Free filter, Amber filter,
Filters			UV Cut (Normal Use), IR Cut (Normal Use)	UV Cut (normal use), IR Cut (normal use),	
		UV Cut (Normal Use), IR Cut (Normal Use)		ND filter (13% transmission)	ND filter (13% transmission), exciter filter for slit lamp*,
				No filed (15% transmission)	barrier filter for slit lamp*
					barrier filter for slit lamp.
Illumination Lamp	LED	SL-D2/D4: Halogen		Halogen	Halogen
		SL-D4 (Type:LED): LED			LED
Base Longitudinal Movement	90mm		90mm	90mm	90mm
Lateral Movement	100mm	100mm		100mm	100mm
Fine Base Movement	12mm	SL-D2/D4: 15mm , SL-D4LED: 12mm		12mm	12mm
Vertical Movement	30mm	30mm		30mm	30mm
vertical Movement	SOMM		3011111	John	Somm
Chinrest					
Vertical Movement	80mm	80mm		80mm	80mm
Dawer					
Power Input (Primary)	100-240V	100-240V		100-240V	100-240V
Frequency	50/60Hz	50/60Hz		50/60Hz	50/60Hz
Power Consumption	40VA	SL-D2/D4:160VA, SL-D4 (Type:LED):110VA		80VA	110VA
Dimensions			1		
w/ Table	550mm(W) x 399mm(D) x 520mm(H)	550mm x 399mm x 558mm		550mm x 430mm x 760-790mm	550mm(W) x 420mm(D) x 752 to 782mm(H)
w/ Unit Table	440mm(W) x 379mm x 558mm(H)	440mm x 379mm x 558mm		440mm x 410mm x 760-790mm	440mm(W) x 397mm(D) x 752 to 782mm(H)
w/o Table and Chinrest	329mm(W) x 293mm(D) x 414-444mm(H)	<b>SL-D2/D4:</b> 329mm x 306mm x 415-445mm,	<b>SL-D4 (Type:LED)</b> 329mm(W) x 331mm(D) x 426 to 456mm(H)	312mm x 296mm x 646-676mm	329mm(W) x 331(304)mm(D) x 652 to 682mm(H)
Weight					
w/ Table	15.2 Kg	17kg		17.5kg	19kg
	-	16kg			401
w/ Unit Table	14.8 Kg		16kg	17kg	18kg

### $^{*}$ A model with standard cobalt blue filter instead of the exciter/barrier filter combination is also available.

### DC-4 Digital Camera Unit (Option)

System	Beam flux division by beam splitter		
Beam splitter division ratio	Operator side: 50%; Camera side: 50%		
Focus distance	35.7mm		
Sensor type	1/2.5 type CMOS; Effective pixel number: 5.04 M pixels		
Dimensions	76mm(W) x 163mm(D) x 304mm(H)		
Weight	430g		



DC-4
Digital Camera Unit



SL-2G SL-D2
Slit Lamp Slit Lamp



SL-D4/D4 (Type:LED)
Slit Lamp



SL-D301 Slit Lamp



SL-D701 Slit Lamp





### M TOPCON CORPORATION

75-1 Hasunuma-cho, Itabashi-ku, Tokyo 174-8580, JAPAN. Phone: +81-(0)3-3558-2522/2502 Fax: +81-(0)3-3965-6898 https://topconhealthcare.jp

### TOPCON SINGAPORE MEDICAL PTE. LTD.

100G Pasir Panjang Road, #05-05, Interlocal Centre, SINGAPORE 118523 Phone:+65-68720606 Fax:+65-67736150

 $\hbox{E-mail: med.sales.sg@topcon.com}$ 

https://topconhealthcare.sg

### TOPCON INSTRUMENTS ( MALAYSIA ) SDN. BHD.

No.6, Jalan Pensyarah U1/28, Hicom Glenmarie Industrial Park, 40150 Shah Alam, Selangor, MALAYSIA Phone: +603 5022 3688 Fax: +603 5031 3968 Email: mys\_medical\_sales@topcon.com

https://topconhealthcare.my

### TOPCON INSTRUMENTS (THAILAND) CO., LTD.

77/162 Sinnsathorn Tower, 37th Floor, Krungthonburi Rd., Klongtonsai, Klongsarn, Bangkok 10600, THAILAND Phone: 662-4401152 Fax: 662-4401158

 $Email: tha\_medical@topcon.com$ https://eyecare.topcon.co.th

### MEHRA EYETECH PRIVATE LIMITED

801 B Wing, Lotus Corporate Park, Graham Firth Steel Compound Goregaon (East) Mumbai 400063 Maharashtra, INDIA Phone: +91-22-61285455 https://topconhealthcare.in

### TOPCON (BEIJING) MEDICAL TECHNOLOGY CO., LTD.

Room 2808, Tower C, JinChang An Building, No.82, Middle Section of East 4th Ring Road, Chaoyang District, Beijing 100124, P.R. CHINA Phone: +86-10-8794-5176 http://www.topcon-china.net/

IMPORTANT In order to obtain the best results with this instrument, please be sure to review all user instructions prior to operation.

©2022 Topcon Corporation

