

Refraction Unit CS-18 Plus



Technical Specification

Motor	: 24 V DC Actuator
Elevation Strokes	: 200 mm
Seat Height	: 505 mm to 705 mm
Chair Up-Down Movement	: With Footswitch
Inclination of Backrest	: No Inclination
Rotation	: No Rotation (Fix Type)
Armrest	: Reclining Type
Footrest	: Yes
Load Lifting Capacity	: 190 Kgs.
Table Top Size	: 405 mm x 635 mm
Table Top Weight Capacity	: 25 Kg. Single Instrument
Table Top Height – From Floor	: 920 mm
Table Top Movement	: 0 to 90 Degree
Base Plate Size	: 655 mm x 650 mm
Trial Lens Tray Drawer	: Sliding Drawer Dust Proof
Light Pole Height From Floor	: 1745 mm
Light Pole Size	: 1093 mm x 42 mm Diameter
Over Head Light	: 5W LED Bulb
AUX Plug on Unit	: 2 Nos. of 220V AC Power Output
Power Source	: 220V AC, 50 HZ, 5 AMP
Maximum Area Required	: 1475 mm (58") x 1100 mm (43.5")
Gross Weight	: 130 Kgs.



58, Tirupati Industrial Estate,
 Nr. Jay Bharat Rangshala Compound,
 B/h. Amber Cinema, Saraspur,
 Ahmedabad-380024 (India)
 Tel. : (079) 2274 8813, 2274 9913
 Fax : (079) 2276 0870.
 E-mail : info@plantechindia.com
 Web : www.plantechindia.com