

Refraction Unit CS-18



Technical Specification

Motor	: 24 V DC Actuator
Elevation Strokes	: 200 mm
Seat Height	: 505 mm to 705 mm
Chair Up-Down Movement	: With Footswitch
Inclination of Backrest	: No Inclination – Fix Type
Rotation	: No Rotation (Fix Type)
Armrest	: Reclining Type
Footrest	: Yes
Load Lifting Capacity	: 190 Kgs.
Table Top Size	: 405 mm x 635 mm
Table Top Weight Capacity	: 25 Kg. Single Instrument
Table Top Height – From Floor	: 920 mm
Table Top Movement	: 0 to 90 Degree
Base Plate Size	: 655 mm x 650 mm
Trial Lens Tray Drawer	: Sliding Drawer Dust Proof
Power Source	: 220V AC, 50 HZ, 5 AMP
Maximum Area Required	: 1475 mm (58") x 1100 mm (43.5")
Gross Weight	: 130 Kgs.