

Pure Sitting Comfort For Smart Doctors



Doctor's Stool

Comfortable,
Convenient,
Competent Consultation

- ▶ An intelligent, maximum-comfort seating support system for the consulting doctor (Shape: circular, Size: 38x8 cm, Material: polyurethane cushion)
- ▶ Smooth, 360-degree rotation.
- ▶ Easy-to-sit, polyurethane seat surface.
- ▶ In-built gas cylinder enables smooth & efficient seat adjustments (range: 52 to 68 cm).
- ▶ Strong yet soft, semispherical, height-adjustable backrest offers optimal elbow, upper body & back support.
- ▶ Stable, epoxy-coated five-pronged base.
- ▶ Sturdy, smooth-surfaced (chromium-plated) bottom ring for feet support.
- ▶ Smooth movements through five high-quality castors.

[Since Plantech continually upgrades its technology and know-how, Design & Specifications might undergo positive changes from time to time. The company reserves the right to do this.]

MANUFACTURED BY
Plantech
MEDICAL SYSTEMS
PRIVATE LIMITED

58, TIRUPATI IND. ESTATE,
NR. JAY BHARAT RANGSHALA COMP.,
B/H. AMBER CINEMA, SARASPUR,
AHMEDABAD - 380 024 (INDIA)
TEL. : (079) 2274 8813, 2274 9913
FAX : (079) 2276 0870
E-mail : pms@icenet.net / info@plantechindia.com
Web : www.plantechindia.com

Comfort. Convenience. Precision. Flexibility. Scientific Co-ordination.