

# A Reliable Techno-Platform For Doctors



The Motorized Instrumentation Table from Plantech is specially engineered for more leg room for both the doctor and the patient, enabling better interactions.

The strong, sturdy base enables mounting of a range of large-sized testing equipment (Perimeter and other Ophthalmic instruments), and its optimal balance ensures uninterrupted efficiency.

The smooth & noiseless, motor-driven, convenient, feathertouch vertical height adjustment is a testimony of the evolved engineering that goes into this table.

The high-quality castors ensure smooth movements and unshakeable stability.

#### TECHNICAL SPECIFICATIONS

Maximum Table Height	870 mm
Minimum Table Height	670 mm
Vertical Adjustment Range	200 mm
Maximum Load	55 kgs
Table Top Size	570 X 400 X 25 mm
Overall Dimensions (Body)	530 X 580 X 870 mm (W) (L) (H)
Weight Nett.	33 kgs.
Power Supply	AC 230V 50/60Hz.
Power Consumption	250 VA

S p a c i o u s . S t r o n g . S t a b l e . S e n s i b l e .

*Customer-Friendly Design;  
Comfort-Oriented Engineering*

## **MOTORIZED INSTRUMENTATION TABLE**



**Plantech**  
MEDICAL SYSTEMS  
PRIVATE LIMITED