



MOTORISED INSTRUMENT TABLE AIT-13 SPECIFICATION

Maximum Height	: 910 mm
Minimum Height	: 710 mm
Vertical Adjustment Range	: 200 mm
Table Top Size	: 1070 mm x 655 mm x 25 mm
Operational Control	: Up/Down Switch
Maximum Area Required	: 1210 mm x 680 mm
Base Movement	: Jerk Free Movement with Castor Locking
Auxiliary Out-Put	: 1 No. of 220V AC Power Out Put
Current	: 230 V A/C – 50 Hz – 5 AMP
Motor	: 24V DC Actuator
Gross Weight	: Kgs

Additional Accessories

- LCD Mount