

# Phoropter Arm



## PHOROPTER ARM SPECIFICATION

<b>Displacement for Up-Down Travelling</b>		: 290 mm
<b>Main Holding Coller ID Size</b>		: Ø 42.2 mm
<b>Pole to Arm Pin Center Distance</b>	<b>Up</b>	: 486 mm
	<b>Middle</b>	: 525 mm
	<b>Down</b>	: 486 mm
<b>Knob Locking</b>		: For Arm Locking & Phoropter Rotation Locking
<b>Gas Syringe Up-Down Working</b>		: Down with Force & Up with Lever Push
<b>Gross Weight</b>		: 6.5 Kgs