


# A Perfect Choice For Vision Testing



## Plantech Vision Testing Drum

- Plantech Vision Testing Drum is well Design Compact and Light Weight.
- Vision Drum can be Mounted on any Refraction Unit OR It can be Mounted on Wall.
- Vision Drum Body is Made from Cold Rolled Mild Steel and Painted with Epoxy Powder Paints.
- Characters of Vision Drum are Designed As Per Formula of Eye Refraction Correction.
- Each Vision Drum having 5 Nos. Charts of Different Languages.
- All Charts are Available Up To 6/4 Correction.
- Vision Drum Having a Colour Test, Dot and Duochrom Test.
- Vision Drum is Operated By Infrared Remote Control & with Cord Remote Control.
- Noise-free Up & Down Movement.
- Over All Dimensions : 375 x 310 x 140 mm
- Gross Weight : 10 Kg.



58, Tirupati Industrial Estate,  
Nr. Jay Bharat Rangshala Compound,  
B/h. Amber Cinema, Saraspur,  
Ahmedabad-380024 (India)  
Tel. : (079) 2274 8813, 2274 9913  
Fax : (079) 2276 0870.  
E-mail : [info@plantechindia.com](mailto:info@plantechindia.com)  
Web : [www.plantechindia.com](http://www.plantechindia.com)