

# A Perfect Partner For Micro Surgery



## Surgeon's Chair SC-16

**Ergonomy**

**Comfort**

**High Precision, Optimum Flexibility**

- Perfect Ergonomics for precision-driven microsurgery.
- Comfortable sitting angles, fit & flexibility.
- Noise-free Up & Down Movement.
- Easy Push Button Foot-switch controlled height adjustment.
- Strong back support with effortless spring-action & Inclination.
- Backrest cushion with easy-to-use height adjustment system.
- Flexible swivel type, yet sturdy arm supports.
- Lockable wheel for Position Fixing.
- Comfortable seat with polyurethane cushion.

### TECHNICAL SPECIFICATIONS

Maximum Height	: 725 mm
Minimum Height	: 525 mm
Vertical Adjustment Range	: 200 mm
Operational Control	: Push Button Foot Switch
Base Movement	: Free Movement on Castor Wheel
Maximum Load	: 150 Kgs.
Backrest	: Continuous Push Back Support
Arm-Rest	: Adjustable by Vertically & Horizontally & Swivel
Arm-Rest Distance	: 510 mm Min / 630 mm Max.
Weight	: 50 Kgs.
Power Supply	: 230 V A/C – 50 Hz